
**kerecis™**

**Omega3 Wound**

## **INSTRUCTIONS FOR USE**

Naturally occurring lyophilized acellular fish skin.

## **INTENDED USE**

Kerecis™ Omega3 Wound maintains and supports a healing environment for the management of wounds, where viable (red) tissue can be seen in the wound base.

The product is used to treat wounds of the following etiologies:

- Diabetic ulcers
- Chronic vascular ulcers
- Venous ulcers
- Trauma wounds (including abrasions, lacerations and skin tears),
- Acute surgical wounds (including debridement, amputation, and donor sites)
- Surgical wounds (dehiscence or failed healing after surgery)
- Imminent failure of split thickness skin graft
- Post-Injection Necrosis

Kerecis™ Omega3 Wound is made from Atlantic Cod, supplied sterile in peel-open packages, and is intended for one-time use.

Not made with rubber latex.

## **CAUTION**

For health care professional use only. These recommendations are designed to serve only as general guidelines. They are not intended to supersede institutional protocols or professional clinical judgment concerning patient care.

## **CONTRAINDICATIONS**

This device is derived from a fish source and should not be used in patients with a known allergy or other sensitivity to fish material.

## **Kerecis™ Omega3 Wound - PRECAUTIONS**

- Remove most of necrotic tissue /biofilm prior to device application.
- Ensure that the product does not overlap the wound edges but only covers the basis of the wound.
- Remove exudate and control bleeding prior to device application.
- Do not reuse or re-sterilize because of cross contamination and loss of physical characteristic of the material respectively.
- Discard unused portions of the device.
- Sterile if package is unopened and undamaged.
- Do not use the device if the package seal has been broken or if handling has caused damage or contamination.
- Do not use in case of known fish allergies.

## **Kerecis™ Omega3 Wound - POTENTIAL COMPLICATIONS**

The following complications are possible. If any of these conditions occur, the device should be removed.

- Infection
- Chronic inflammation (Initial application of wound dressings may

be associated with transient, mild, localized inflammation)

- Allergic reaction
- Excessive redness, pain, swelling, or blistering

## **STORAGE**

This device should be stored in a clean, dry location at room temperature.

## **STERILIZATION**

This device has been sterilized using ethylene oxide sterilization.

## **SUGGESTED INSTRUCTIONS FOR USING Kerecis™ Omega3 Wound**

The following instructions are designed to serve only as general guidelines. They are not intended to supersede institutional protocols or professional clinical judgment concerning patient care. Always handle Kerecis™ Omega3 Wound in an aseptic manner.

## **INITIAL APPLICATION**

1. When applying Kerecis™ Omega3 Wound for the first time, clean the wound bed by removing necrotic tissue/biofilm to obtain a fresh tissue surface. Control bleeding.
2. Remove the Kerecis™ Omega3 Wound sheet from the pouch in an aseptic manner.
3. Cut the sheet roughly to the size of the wound so that the sheet will not cover the wound edges and pre-hydrate in sterile saline (0,9% w/v sodium chloride Ph.Eur) for 30-180 seconds.
4. Apply the sheet directly into the wound with its outer side, identifiable by its scaly-texture, facing upwards. Make sure that the

sheet does not overlap the wound edges.

5. More than one sheet may be necessary for complete coverage of large wounds. Overlap wound sheet edges slightly to assure coverage of the entire wound.
6. Apply an appropriate secondary non-adherent wound dressing to maintain a moist wound environment. Be careful to use a fixation that does not soak the wound, e.g., a non-breathable plastic film.

### **FOLLOW-UP APPLICATION**

After starting wound management with Kerecis™ Omega3 Wound the wound needs to be checked regularly to see if the device needs to be reapplied, and if the secondary dressing is maintaining a sufficiently moist wound environment.

1. Inspect the wound initially every two days. The duration between inspections can be extended up to seven days as healing progresses.
2. Insert a new Kerecis™ Omega3 Wound sheet if the previously applied sheet has been absorbed and is no longer visible.
3. The device is intended to stay in the wound bed. If removal of remaining Kerecis™ Omega3 Wound is deemed desirable it should be gently lifted with tweezers out of the wound bed.
4. Change the secondary dressing as needed to maintain a moist, clean wound area.

As wound healing occurs, redness and swelling of wound edges will decrease and the level of exudate will be reduced. These are signs of wound healing and are often seen before epithelization is obvious.

Always make sure to use an appropriate type of secondary wound dressing to maintain a moist, clean wound area.

Should you suspect any adverse reaction from this product please report it to the Email [adversereactions@kerecis.com](mailto:adversereactions@kerecis.com) leaving a phone number and Email address. Adverse reactions might also be subject to mandatory reporting to the authorities. Please submit the data on the package or pouch (e.g., "lot number" and "use before" date).

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Kerecis Limited, Eyrargata 2, 400 Isafjordur, Iceland. Telephone: +354 562 2601, Email: [info@kerecis.com](mailto:info@kerecis.com), [www.kerecis.com](http://www.kerecis.com). Patented: US 8.613.957. Kerecis is a trademark or registered trademark of Kerecis Limited.