Kerecis

Kerecis
MariSurgi™

EN - INSTRUCTIONS FOR USE

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INSTRUCTIONS FOR USE

Fish skin for Surgical Use

Kerecis MariSurgi™

MariSurgi™ is an acellular fish skin medical device for wound management available in a variety of sizes and configurations. It is made of decellularized and lyophilized fish skin of the North Atlantic cod (*Gadus morhua*) and supplied sterile in Tyvek® peelable pouch packaging.

MariSurgi™ supports the body's own abilities to regenerate damaged tissue. MariSurgi™ aids the efficient closure of acute and chronic wounds through earlier development of granulation tissue and healing by secondary intent.

INTENDED USE AND INDICATIONS

MariSurgi™ can be used on patients, 18 years or older for a variety of wound types, including:

- Partial and full-thickness wounds, including:
 - Diabetic foot ulcers
 - Venous ulcers
 - Arterial ulcers
 - Partial-thickness burns
- Trauma wounds (including abrasions, lacerations, skin tears)
- Acute and chronic surgical wounds (including post-amputation, dehiscence, and auto-graft donor sites)

CAUTION

For use by a licensed healthcare professional, only trained in wound management and care.

CONTRAINDICATIONS

Do not use MariSurgi™ in patients with known history of allergy or hypersensitivity to fish or fish products.

Shellfish allergy is not a contraindication.

Do not apply MariSurgi™ over infected or deteriorating wounds until the underlying cause has been resolved.

WARNING AND PRECAUTIONS

Single patient use only; do not reuse or re-sterilize the device.

- Sterile only if the package is unopened and

undamaged.

- Do not use the device if the package seal has been broken, unintentionally opened, or if handling caused damage or contamination. Dispose of such device following institutional protocols.
- Do not use past the expiration date as indicated on the product label.
- MariSurgi™ is not designed, sold, or intended for use except as indicated.
- No clinical trials have been conducted in pediatric, pregnant and breastfeeding populations.

RESIDUAL RISK

See Sections: Warnings and Precautions, Contraindication, and Potential Complications.

INSTRUCTION FOR USE

The instructions serve as general guidelines, only. Applicable to solid and meshed variations of the device. The instruction for use should not supersede institutional protocols or professional judgment concerning patient care.

MariSurgi™ should always be handled using aseptic technique and applied onto viable tissue, only.

First – Double Pouch Handling, *where applicable*

1. Refer to label if device is supplied in a single or double pouch sterile barrier (see Symbols from Device Packaging).
2. Open outer pouch within appropriate environment.

3. Transfer inner pouch to sterile field using aseptic technique.

Second – Wound Preparation

Prior to applying the product, prepare the wound bed using standard debridement techniques. Establish a wound free of debris, exudate, and non-viable tissue. Irrigate and control significant bleeding. Post-debridement pinpoint bleeding often indicates adequate debridement and may be beneficial in ensuring optimal take.

Third – Initial Application

1. Remove MariSurgi™ from the pouch in an aseptic manner.
2. If needed cut MariSurgi™ to fit the wound size with minimal overlap at the wound edge.
3. Hydrate the device in room-temperature sterile saline for at least 30 seconds but not longer than 10 minutes for improved pliability during application. Further trimming of the hydrated device is possible, if required.

It is recommended to stretch the meshed device variant at least once before application to the wound bed.

4. Place the device directly into the wound bed. When using multiple sheets, overlap sheet edges slightly to ensure full coverage.
5. Secure MariSurgi™ using for example adhesive plaster strips, a silicone-based contact layer, surgical staples, sutures, or skin glue. Bolster or use VAC to ensure firm contact of the product with the wound bed.
6. Dispose of any unused material of the device following institutional protocols.

Fourth – Final Step

Use appropriate cover/secondary dressing(s) to maintain a moist wound environment.

Ensure appropriate exudate management. Change cover/secondary dressing(s) as needed.

The recommended secondary dressings include paraffin gauze (ex. Jelonet®), semi-permeable films (ex. 3M™ Tegaderm™ Transparent Film Dressing), foams (ex. Allevyn Foam, UrgoStart® Dressings, TheraSorb®), hydrofiber dressings (ex. AQUACEL® EXTRA™ Dressing), hydrogels (HydroClean®), alginates (Algostéril) and cellulose based dressings (ex. 3M™ Adaptic™ Non-Adhering Dressing).

Avoid non-permeable dressings to prevent maceration of the wound.

Follow-Up and Subsequent Applications

Ensure an optimal wound healing environment by regular monitoring of the wound after the initial application of MariSurgi™.

- Inspect wound 2 to 3 days after initial application. The duration between inspections can be extended up to 7 days as healing progresses. Change the cover/secondary dressing(s) as needed.
- The medical device will degrade in the wound. Additional applications of the device can be performed if the previously applied device is no longer visible in the wound.

POTENTIAL COMPLICATIONS POST APPLICATION OF THE DEVICE

Monitor the wound area for any of the following:

- Infection
- Chronic inflammation
- Allergic reaction
- Moderate to severe redness (when compared to initial application), pain, swelling or blistering

If any of the above mentioned occur the device should be removed by a licensed healthcare professional in an aseptic manner.

STORAGE

Store at room temperature in a clean and dry location away from direct sunlight.

STERILIZATION

This device has been sterilized using ethylene oxide.

PACKAGING

MariSurgi™ is supplied sterile in single or double pouch peel-open packaging.

DISPOSAL

The device should be disposed according to institutional protocols for animal-derived tissue.

ADVERSE REACTION / SERIOUS INCIDENT

Any adverse reaction or a serious incident suspected from this device should be reported to adversereactions@kerecis.com. Please include at least a contact phone number, Email address, and the device lot number and expiration date. Note that adverse reactions and serious incidents are subject to mandatory reporting to the competent authorities having jurisdiction.

LINK ON EUDAMED UDI DATABASE

SSCP is available in the European database on medical devices (Eudamed), linked to the basic UDI-DI

- URL to the Eudamed public website: <https://ec.europa.eu/tools/eudamed>
- Basic UDI-DI 569431096FishSkinNG

SYMBOLS FROM DEVICE PACKAGING

	MEDICAL DEVICE		DISTRIBUTOR
	CONTAINS BIOLOGICAL MATERIAL OF ANIMAL ORIGIN		DO NOT RE-STERILIZE
	KEEP AWAY FROM SUNLIGHT		SINGLE STERILE BARRIER SYSTEM
	KEEP DRY		DOUBLE STERILE BARRIER SYSTEM
	UPPER LIMIT OF TEMPERATURE, 75°F, 25°C, LOWER LIMIT OF TEMPERATURE, 32°F, 0°C		MANUFACTURER
	CONSULT INSTRUCTION FOR USE		CATALOG NUMBER
	DO NOT USE IF PACKAGING IS DAMAGED		UNIQUE DEVICE IDENTIFIER
	DO NOT REUSE		LOT NUMBER
	CAUTION		EXPIRY DATE
	STERILIZED USING ETHYLENE OXIDE		CE MARK, NOTIFIED BODY: DEKRA CERTIFICATION B.V.

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trademarks of Coloplast A/S
MARCADO CE, CUERPO NOTIFICADO: DEKRA CERTIFICATION B.V.
Notified Body: DEKRA Certification B.V.