


GraftGuide®

Fish skin application for burns




GraftGuide®

GraftGuide is intact fish skin used to regenerate tissue in surgical, traumatic, and acute wounds.³ Because there is no known risk of viral or prion disease transfer from North Atlantic cod to humans, the fish skin only needs mild processing with our proprietary method.²

This mild processing preserves the skin's natural qualities², including its three-dimensional structure, mechanical properties, molecular organization, and composition.¹⁻⁴

GraftGuide promotes healing with minimal impairment of functionality and positive cosmetic outcomes.⁴ The product is homologous to human skin¹ and when applied to damaged tissue, such as a burn or wound, helps to support the body's own cells to regenerate tissue.⁵⁻¹²



GraftGuide Meshed 3-to-1 is designed to easily expand across large, partial-thickness burns while preventing the rebound effect sometimes experienced when grafts with lesser expansion ratios are stretched across large burns.

This reduced rebound effect combined with strength and durability makes affixation easier and quicker for large surface area wounds. The larger expansion ratio also helps reduce "tenting" that can occur when 2-to-1 pre-meshed grafts or solid grafts are applied to burns with significantly uneven wound beds.



GraftGuide® Micro



GraftGuide Micro is intact fish skin that has been fragmented. It offers more surface area than the non-fragmented version, and is designed to adhere to, and fill, deep wound spaces and irregular geometries.



Intended Use

The product is indicated for the management of wounds including:

- Partial and full-thickness wounds
- Trauma wounds including second-degree burns
- Surgical wounds

Features and Benefits

GraftGuide pioneers the use of fish skin for the management of burn wounds.

- Provides a natural microbial barrier¹, protecting against environmental insults
- Preserves the structure and porosity of the fish skin graft¹, ideal for cellular ingrowth^{1,10}
- Adheres to burn wounds in a variety of shapes and sizes, including mobile joints
- Resists shearing forces, allowing early physiotherapy¹⁰
- Supports the restoration of the skin's pigmentation, while minimizing scarring^{4,6,10-12,15}

GraftGuide® Mano



GraftGuide Mano is intact fish skin specifically designed to reduce operation time when managing burn injuries, with the product's unique shape allowing it to easily cover the complex 3D structure of the hand. It is available in two sizes, medium and large, and in either left or right versions, for use in palmar or dorsum applications.



Kerecis Mechanism of Action

Kerecis patented fish skin is only gently processed, preserving its homologous structure to human skin. This product has had its cells removed and a native structure largely unaffected by processing, both traits recognized as optimal for a tissue replacement product.¹³

INTACT MOLECULAR ORGANIZATION
Fish skin is homologous to human skin²


THREE-DIMENSIONAL STRUCTURE
Intact structure provides a framework to support tissue regeneration with limiting scarring¹⁴

NATURAL MECHANICAL PROPERTIES
Naturally strong, handles like skin, and easy to suture or staple³

PRESERVED MOLECULAR CONTENT
Chemical complexity of the fish skin supports rapid cell ingrowth and neovascularization¹


Homologous to Human Skin

FISH SKIN



100 μm

HUMAN SKIN



100 μm


Epidermis

Dermis

Largest RCT on DFUs with Exposed Bone or Tendon Brings New Hope for Patients

255 patients, 15 centers, UT Grades 2 & 3


Proportion of Wound Closure: IFSG vs. SOC



Intact fish skin heals 66% more wounds than SOC at 16 weeks

Intact Fish Skin Graft (Orange)
Standard of Care (Dark Blue)

The Use of Intact Fish Skin as a Novel Treatment Method for Deep Dermal Burns following Enzymatic Debridement: A retrospective case-control study




A) There was a significant reduction in wound size over time in fish skin treated wounds compared to Suprathel-treated wounds ($p < 0.001$). $N = 12$.
B) The application of the definitive wound closure to the point of 95% epithelialization: Suprathel (45.6 ± 6.6 days), STSG (34.7 ± 12.5 days), and fish skin graft (22 ± 6.3 days).

CONCLUSIONS

Enzymatic debridement in combination with intact fish skin grafts resulted in improved short term and long term outcomes.

- Faster burn healing
- Reduced need for autografting
- Better functional outcomes
- Better scar outcomes



GraftGuide Cases

Healing by secondary intent with epithelialization of the skin resulting in desirable cosmetic and functional outcomes.^{7,15}


Case courtesy of Dr. Alfredo Cordova, MD, FACS, FAAST

93% TBSA Full Thickness Flame Wound

Background: 28 y/o Male with Schizophrenia. Poured gasoline and set himself on fire.

Previous Management: Allograft

Applications: Kerecis Meshed 2:1 Graft



Case courtesy of Jason M. Bregg, MD

Deep Partial-Thickness Hand Burn

Previous Management: Presented to the burn center on the same day as the injury with no previous treatment reported.

Applications: GraftGuide Mano (Size L)




Case courtesy of Ram Velamuri, MD, FACS, MS, DNB, MRCS-Ed

40% TBSA Burn

Previous Management: BTM was applied but removed due to infection that came up.

Applications: Graftguide Meshed 2:1



GraftGuide®

SOLID

Catalog # Box of 10	Catalog # Single Unit	Product Size	Coverage
50202S03D2D	50202S03D0D	7 x 10 cm	70
50202S21D2D	50202S21D0D	7 x 20 cm	140
50202S24D2D	50202S24D0D	250 cm ²	250
50202S23D2D	50202S23D0D	300 cm ²	300

MESHED 2:1

Catalog # Box of 10	Catalog # Single Unit	Product Size	Coverage
50202N03D2D	50202N03D0D	7 x 10 cm meshed 2:1	up to 126
50202N21D2D	50202N21D0D	7 x 20 cm meshed 2:1	up to 252
50202N24D2D	50202N24D0D	250 cm ² meshed 2:1	up to 450
50202N23D2D	50202N23D0D	300 cm ² meshed 2:1	up to 540

MESHED 3:1

Catalog # Box of 10	Catalog # Single Unit	Product Size	Coverage
50202023B2D	50202023B0D	300 cm ² meshed 3:1	up to 840

GraftGuide® Micro

Catalog # Box of 10	Catalog # Single Unit	Product Size	Coverage
50202P32D2D	50202P32D0D	114 cm ²	up to 114

GraftGuide® Mano

Catalog # Box of 10	Catalog # Single Unit	Product Type / Size	Coverage
50202N33D0D	50202N33D2D	Dorsum L/Palmar R / M	147
50202N34D0D	50202N34D2D	Dorsum R/Palmar L / M	147
50202N35D0D	50202N35D2D	Dorsum L/Palmar R / L	196
50202N36D0D	50202N36D2D	Dorsum R/Palmar L / L	196

Sustainably Harvested in Iceland


Kerecis is pioneering the use of fish skin in the globally expanding cellular therapy and regenerative medicine market. It is the only regulatory cleared manufacturer of fish-skin-based medical devices for tissue regeneration in the world.

Kerecis fish skin is a byproduct of fish processing in Iceland. The fish are caught from sustainable wild fishing stocks in the North Atlantic Ocean.

The Kerecis ISO 13485 manufacturing facilities are in the town of Isafjordur in the Westfjords region, 20 miles south of the Arctic Circle, and are powered by renewable energy.




KERECIS PRODUCT CONFIGURATIONS



Solid




Fenestrated




Fragmented



Meshed 2:1



Meshed 3:1



Integral Silicone

REFERENCES

- 1) Magnusson, S., et al. Decellularized fish skin: characteristics that support tissue repair. *Laeknabladid*. 101, 567-573 (2015).
- 2) Rakers, S., et al. 'Fish matters': the relevance of fish skin biology to investigative dermatology. *Exp Dermatol*. 19, 313-324 (2010).
- 3) Baldursson, B. T., et al. Healing rate and autoimmune safety of full-thickness wounds treated with fish skin acellular dermal matrix versus porcine small-intestine submucosa: a noninferiority study. *Int J Low Extrem Wounds*. 14, (2015).
- 4) Kirsner, R. S., et al. Fish skin grafts compared to human amnion/chorion membrane allografts: A double-blind, prospective, randomized clinical trial of acute wound healing. *Wound Repair and Regeneration*, 28(1), 75-80 (2020).
- 5) Lullove, E. J., et al. A multicenter, blinded, randomized controlled clinical trial evaluating the effect of Omega-3-rich fish skin in the treatment of chronic, nonresponsive diabetic foot ulcers. *Wounds* (2021).
- 6) Winters, C., et al. Cost Effectiveness of Fish Skin Grafts Versus Standard of Care on Wound Healing of Chronic Diabetic Foot Ulcers: A Retrospective Comparative Cohort Study. *Wounds*. 32(10):283-290.(2020).
- 7) Stone, II, R., et al. Accelerated Wound Closure of Deep Partial Thickness Burns with Acellular Fish Skin Graft. *Int J Mol Sci*. 22(4):1590 (2021).
- 8) Yang, C.K., et al. A Prospective, Postmarket, Compassionate Clinical Evaluation of a Novel Acellular Fish-skin Graft Which Contains Omega-3 Fatty Acids for the Closure of Hard-to-heal Lower Extremity Chronic Ulcers. *Wounds*. 28(4):112-8 (2016).
- 9) Trinh, T. T., et al. Marine Omega3 Wound Matrix for the Treatment of Complicated Wounds. *Phlebologie*, 45, 93-98 (2016).
- 10) Woodrow, T., et al. Treatment of diabetic foot wounds with acellular fish skin graft rich in omega-3: a prospective evaluation. *J Wound Care*. 28, 76-80 (2019).
- 11) Sitje, T. S., et al. A Clinical innovation: fish-derived wound product for cutaneous wounds. *Wounds Int*. 2018 9, 44-50 (2018).
- 12) Patel, M. & Lantis II, J. C. Fish skin acellular dermal matrix: potential in the treatment of chronic wounds. *Chronic Wound Care Manag Res*. 6, 59-70 (2019).
- 13) Sibbald, R., et al. Special considerations in wound bed preparation 2011: An update. *Adv Skin Wound Care*. 24(9):415-36 (2011).
- 14) Wallner, C., et al. The Use of Intact Fish Skin as a Novel Treatment Method for Deep Dermal Burns Following Enzymatic Debridement: A Retrospective Case-Control Study. *Euro Burn J*. 3(1):43-55 (2022).
- 15) Andrews JW, Adams MJ, Montenegro-Johnson TD. A universal scaling law of mammalian touch. *Sci Adv*. 2020 Oct 9;6(41):eabb6912. doi: 10.1126/sciadv.abb6912. PMID: 33036967; PMCID: PMC7546702.

kerecis[®]

OUR VISION
To extend life by supporting
the body's own ability
to regenerate

U.S. and international
patents and trademarks
granted and pending

KERECIS UNITED STATES
2101 Wilson Blvd., Suite 900
Arlington, Virginia 22201
+1 (703) 287-8752
vet@kerecis.com

www.kerecis.com