

Burn Cases

Fish skin application for burn wounds



About Kerecis

Kerecis fish skin is homologous to human skin¹⁷, and is used to support tissue regeneration and repair.² Kerecis fish skin products are FDA cleared for multiple clinical applications.

Because no known viral transfer risk exists between Atlantic Cod and humans, the fish skin needs only mild processing. Our proprietary manufacturing method preserves the skin's natural qualities, including its three-dimensional structure, mechanical properties, molecular organization, and composition.^{2,3}

When grafted onto damaged human tissue, such as a burn or a wound, the fish skin recruits the body's own cells, supporting the ability to regenerate.¹

The superior clinical and economic performance of Kerecis fish skin has been demonstrated in multiple blinded, randomized, controlled clinical trials^{3,4,5} and numerous other clinical studies.^{6-13,15} Since there are no known religious or cultural barriers associated with Kerecis products, they can help diverse communities.¹⁴

CONTRAINDICATIONS

Kerecis should not be used in patients with known fish allergies.

Kerecis Mechanism of Action

Key factors for an optimal tissue replacement product are recognized to be a product that has had its cells removed, and a native structure that is not affected or altered by the applied treatment process, thereby retaining its overall lipid content, ECM biochemical composition, structure, and complexity.¹



THREE-DIMENSIONAL STRUCTURE

Intact structure provides a framework to support tissue regeneration¹⁶



NATURAL MECHANICAL PROPERTIES

Naturally strong, handles like skin, and easy to suture or staple³



PRESERVED MOLECULAR CONTENT


Chemical complexity of the fish skin provides the ideal environment for wound healing¹



INTACT MOLECULAR ORGANIZATION

Fish skin is homologous to human skin¹⁷


Fish Skin



Epidermis

Dermis

Human Skin



DISCLAIMER

This communication is to serve only as general information or guidance. Refer to packet inserts for Indications, Contraindications, Cautions and Instructions. None are intended to supersede institutional protocols or professional clinical judgment concerning patient care.

CASE FROM

Dr. Alfredo Cordova, MD, FACS, FFAST

PATIENT PRESENTATION

- 28-year-old male with schizophrenia
- Poured gasoline on himself and set himself on fire
- Sustained 93% total body surface area (TBSA) full-thickness burns

PREVIOUS MANAGEMENT

Allograft

APPLICATION

Kerecis Meshed 2:1 Graft

CLINICAL OUTCOME

Intact fish skin integration and adequate granulation tissue was evidenced in 95% of the surface area as early as 10 days after product application. Patient was later resurfaced with widely meshed (6:1) split thickness autograft and cultured epithelial autograft (CEA).




INITIAL PRESENTATION 93% TBSA burns sustained. Needed to be resuscitated



DAY 0 Underwent fascial excision and initially resurfacing with allograft to his anterior trunk and bilateral lower extremities. Concurrent Management™ BTM, RECELL, Meeks micrograting 9:1, Suprathel



DAY 0 Allografts became slimy and heavily colonized



DAY 0 Allografts removed and Kerecis was placed



DAY 15 Multiple Kerecis 300cm² Meshed 2:1 grafts applied



DAY 27 Good epithelialization. Wound beds optimal for grafting, widely meshed (6:1) split thickness autograft and cultured epithelial autograft applied and healing with no complications

Scan QR-Code for case presentation



CASE FROM

Jason M. Bregg, MD

PATIENT PRESENTATION

- 39-year-old male
- Presented with 4% total body surface area (TBSA) deep partial-thickness thermal burn

PREVIOUS MANAGEMENT

Presented to the burn center on the same day as the injury with no previous treatment reported.

APPLICATION

The wound was debrided using dermabrasion, and GraftGuide Mano (Size L) was applied to cover the dorsal aspect of the right hand. The product was secured in place, and a secondary dressing was applied.

CLINICAL OUTCOME

Single application of GraftGuide Mano resulted in good secondary intention healing for the patient without the need of autograft. Good functional outcome with no evidence of scar contracture.



INITIAL PRESENTATION Same day, post burn injury, prior to fish skin application



DAY 2
Pre-debridement



DAY 2
Excisional debridement with dermatome



INTRA-OP
Application of GraftGuide Mano



DAY 5 Product has almost incorporated fully into the wound bed. No signs of infection



DAY 17 Full epithelialization without the need for a split thickness skin graft

CASE FROM

Ariel Aballay, MD

PATIENT PRESENTATION


- 35-year-old male
- Previous medical history of Afib
- Patient experienced a high voltage electrical injury (7200 V) at work.
- He experienced 15% TBSA burns, mostly full thickness with severe loss of tissue on the left shoulder

APPLICATION

Kerecis Meshed 2:1 Graft was applied to the wound. A biopsy was taken and the wounds were redressed.

CLINICAL OUTCOME

The use of fish skin allowed this patient to progress to healing. Note that provider performed a biopsy which indicated structure formed within fish skin graft is indistinguishable from naturally occurring granulation tissue.



INITIAL PRESENTATION



Kerecis Meshed 2:1 Graft Application



5 days post application



FINAL PRESENTATION

CASE FROM

Richard Ehrlichman, MD

PATIENT PRESENTATION

- 67-year-old male
- Patient sustained a dorsal foot burn 3 weeks before admission presented with cellulitis
- The 13 x 12 cm wound was excised and covered with Kerecis meshed 2:1 Graft

APPLICATION

Kerecis Meshed 2:1 Graft was placed over the wound with a VAC.

CLINICAL OUTCOME

The wounds covered by Kerecis Meshed 2:1 Graft remained clean, the patient was comfortable.

There was excellent granulation tissue presented without infection.



INITIAL PRESENTATION



DAY 1

Application of Kerecis Meshed 2:1 Graft



5 days post application



FINAL PRESENTATION

16 days later

CASE FROM

Ram Velamuri, MD, FACS, MS, DNB, MRCS-Ed

PATIENT PRESENTATION

- 51-year-old male
- 3rd degree burn
- 40% TBSA

APPLICATION

Graftguide Meshed 2:1 was applied

CLINICAL OUTCOME

Patient was grafted and had a small pressure ulcer due to the ACE wraps that were used in rehab. One of his limbs was salvaged and the other was amputated but is able to successfully use his prosthetic and walk.



INITIAL PRESENTATION 40% TBSA



WEEK 1 BTM application. Patient developed infection and BTM product was removed. Kerecis was applied following this image.



2 WEEKS POST FIRST APPLICATION
Granulation tissue developed



FINAL PRESENTATION
1 month post STSG

Scan QR-Code for
case presentation



Sustainably Harvested in Iceland


Kerecis is pioneering the use of fish skin in the globally expanding cellular therapy and regenerative medicine market. It is the only approved manufacturer of fish-skin-based medical devices for tissue regeneration in the world.

Kerecis fish skin is a byproduct of fish processing in Iceland. The fish are caught from sustainable wild fishing stocks in the North Atlantic Ocean.

The Kerecis ISO 13485 manufacturing facilities are in the town of Isafjordur in the Westfjords region, 20 miles south of the Arctic Circle, and are powered exclusively by renewable energy.



KERECIS PRODUCT VARIATIONS




Solid




Fenestrated



Fragmented



Meshed 2:1



Integral Silicone

REFERENCES

- Magnusson, S., et al. Decellularized fish skin: characteristics that support tissue repair. *Laeknabladid*. 101, 567–573 (2015).
- Rakers, S., et al. 'Fish matters': the relevance of fish skin biology to investigative dermatology. *Exp Dermatol*. 19, 313–324 (2010).
- Baldursson, B. T., et al. Healing rate and autoimmune safety of full-thickness wounds treated with fish skin acellular dermal matrix versus porcine small-intestine submucosa: a noninferiority study. *Int J Low Extrem. Wounds*. 14, (2015).
- Kirsner, R. S., et al. Fish skin grafts compared to human amnion/chorion membrane allografts: A double-blind, prospective, randomized clinical trial of acute wound healing. *Wound Repair and Regeneration*, 28(1), 75–80 (2020).
- Lullove, E. J., et al. A multicenter, blinded, randomized controlled clinical trial evaluating the effect of Omega-3-rich fish skin in the treatment of chronic, nonresponsive diabetic foot ulcers. *Wounds* (2021).
- Winters, C., et al. Cost Effectiveness of Fish Skin Grafts Versus Standard of Care on Wound Healing of Chronic Diabetic Foot Ulcers: A Retrospective Comparative Cohort Study. *Wounds*. 32(10):283-290.(2020).
- Stone, Il, R., et al. Accelerated Wound Closure of Deep Partial Thickness Burns with Acellular Fish Skin Graft. *Int J Mol Sci*. 22(4):1590 (2021).
- Yang, C.K., et al. A Prospective, Postmarket, Compassionate Clinical Evaluation of a Novel Acellular Fish-skin Graft Which Contains Omega-3 Fatty Acids for the Closure of Hard-to-heal Lower Extremity Chronic Ulcers. *Wounds*. 28(4):112-8 (2016).
- Trinh, T. T., et al. Marine Omega3 Wound Matrix for the Treatment of Complicated Wounds. *Phlebologie*. 45, 93–98 (2016).
- Woodrow, T., et al. Treatment of diabetic foot wounds with acellular fish skin graft rich in omega-3: a prospective evaluation. *J Wound Care*. 28, 76–80 (2019).
- Sitje, T. S., et al. A Clinical innovation: fish-derived wound product for cutaneous wounds. *Wounds Int*. 2018 9, 44–50 (2018).
- Patel, M. & Lantis Il, J. C. Fish skin acellular dermal matrix: potential in the treatment of chronic wounds. *Chronic Wound Care Manag Res*. 6, 59–70 (2019).
- Sibbald, R., et al. Special considerations in wound bed preparation 2011: An update. *Adv Skin Wound Care*. 24(9):415-36 (2011).
- Easterbrook, C. & Maddern, G. Porcine and bovine surgical products: Jewish, Muslim, and Hindu perspectives. *Arch Surg*. 143(4), 366–370 (2008).
- Wallner, C., et al. The Use of Intact Fish Skin as a Novel Treatment Method for Deep Dermal Burns Following Enzymatic Debridement: A Retrospective Case-Control Study. *Euro Burn J*. 3(1):43-55 (2022).
- Wang, Y. et al. Burn injury: challenges and advances in burn wound healing, infection, pain and scarring. *Advanced drug delivery reviews*, 123, 3–17 (2018)
- Magnusson, S. et al. Regenerative and Antibacterial Properties of Acellular Fish Skin Grafts and Human Amnion/Chorion Membrane: Implications for Tissue Preservation in Combat Casualty Care. *Mil.Med*. 182, 383–388 (2017).

kerecis[®]

OUR VISION
To extend life by supporting the body's own ability to regenerate

U.S. and international patents and trademarks granted and pending.

KERECIS UNITED STATES
2101 Wilson Blvd
Suite 900
Arlington, Virginia 22201
Phone: (703) 287-8752
Email: info@kerecis.com

www.kerecis.com