## INSTRUCTIONS FOR USE

Fish skin for Surgical Use

### Kerecis MariSurgi™

MariSurgi™ is an acellular fish skin medical device for wound management available in a variety of sizes and configurations. It is made of decellularized and lyophilized fish skin of the North Atlantic cod (*Gadus morhua*) and supplied sterile in Tyvek® peelable pouch packaging.

MariSurgi™ supports the body's own abilities to regenerate damaged tissue. MariSurgi™ aids the efficient closure of acute and chronic wounds through earlier development of granulation tissue and healing by secondary intent.

## INTENDED USE AND INDICATIONS

MariSurgi™ is indicated for wound types, including:

- Partial and full-thickness wounds, including:
  - Diabetic foot ulcers
  - Venous ulcers
  - Arterial ulcers
  - Pressure ulcers
  - Partial-thickness burns
- Trauma wounds (including abrasions, lacerations, skin tears)
- Acute and chronic surgical wounds (including post-amputation, dehiscence, and auto-graft donor sites)

## CAUTION

For use by a licensed healthcare professional, only trained in wound management and care.

## CONTRAINDICATIONS

Do not use MariSurgi™ in patients with known history of allergy or hypersensitivity to fish or fish products.

Shellfish allergy is not a contraindication.

Do not apply MariSurgi™ over infected or deteriorating wounds until the underlying cause has been resolved.

## WARNINGS AND PRECAUTIONS

Single patient use only; do not reuse or re-sterilize the device.

- Sterile only if the package is unopened and undamaged.
- Do not use the device if the package seal has been broken, unintentionally opened, or if handling caused damage or contamination. Dispose of such device following institutional protocols.
- Do not use past the expiration date as indicated on the product label.
- MariSurgi™ is not designed, sold, or intended for use except as indicated.
- No clinical trials have been conducted in pediatric, pregnant and breastfeeding populations.
- MariSurgi™ must always be covered with an appropriate secondary cover dressing following application to the wound. The secondary dressing should be selected and managed based on wound characteristics and clinical judgment.

## RESIDUAL RISK

See Sections: Warnings and Precautions, Contraindication, and Potential Complications.

## INSTRUCTION FOR USE

The instructions serve as general guidelines, only. Applicable to solid, fenestrated, fragmented and meshed variations of the device. The instruction for use should not supersede institutional protocols or professional judgment concerning patient care.

MariSurgi™ should always be handled using aseptic technique and applied onto viable tissue, only.

### First – Double Pouch Handling, where applicable

1. Refer to label if device is supplied in a single or double pouch sterile barrier (see Symbols from Device Packaging).
2. Open outer pouch within appropriate environment.
3. Transfer inner pouch to sterile field using aseptic technique.

### Second – Wound Preparation

Prior to applying the product, prepare the wound bed using standard debridement techniques. Establish a wound free of debris, exudate, and non-viable tissue. Irrigate and control significant bleeding. Post-debridement pinpoint bleeding often indicates adequate debridement and may be beneficial in ensuring optimal take.

### Third – Initial Application

1. Remove MariSurgi™ from the pouch in an aseptic manner.
2. If needed cut MariSurgi™ to fit the wound size with minimal overlap at the wound edge.
3. Hydrate the device in room-temperature sterile saline for at least 30 seconds but not longer than 10 minutes for improved pliability during application. Further trimming of the hydrated device is possible, if required.

***It is recommended to stretch the meshed device variant at least once before application to the wound bed.***

4. Place the device directly into the wound bed. If using multiple sheets, overlap sheet edges slightly to ensure full coverage.
5. Secure MariSurgi™ using for example adhesive plaster strips, a silicone-based contact layer, surgical staples, sutures, or skin glue. Bolster or use VAC to ensure firm contact of the product with the wound bed.
6. Dispose of any unused material of the device following institutional protocols.

### Fourth – Final Step

Use appropriate cover/secondary dressing(s) to maintain a moist wound environment.

Ensure appropriate exudate management. Change cover/secondary dressing(s) as needed.

Avoid non-permeable dressings to prevent maceration of the wound.

## Follow-Up and Subsequent Applications

Ensure an optimal wound healing environment by regular monitoring of the wound after the initial application of MariSurgi™.

- Inspect wound 2 to 3 days after initial application. The duration between inspections can be extended up to 7 days as healing progresses. Change the cover/secondary dressing(s) as needed.
- The medical device will degrade in the wound. Additional applications of the device can be performed if the previously applied device is no longer visible in the wound.

## POTENTIAL COMPLICATIONS POST APPLICATION OF THE DEVICE

Monitor the wound area for any of the following:

- Infection
- Chronic inflammation
- Allergic reaction
- Moderate to severe redness (when compared to initial application), pain or swelling

If any of the above mentioned occur the device should be removed by a licensed healthcare professional in an aseptic manner.

## STORAGE

Store at room temperature in a clean and dry location away from direct sunlight.

## STERILIZATION

This device has been sterilized using ethylene oxide.

## PACKAGING

MariSurgi™ is supplied sterile in single or double pouch peel-open packaging.

## DISPOSAL

The device should be disposed of according to institutional protocols for animal-derived tissue.

## ADVERSE REACTION / SERIOUS INCIDENT

Any adverse reaction or a serious incident suspected from this device should be reported to [adversereactions@kerecis.com](mailto:adversereactions@kerecis.com). Please include at least a contact phone number, Email address, and the device lot number and expiration date. Note that adverse reactions and serious incidents are subject to mandatory reporting to the competent authorities having jurisdiction.

## SYMBOLS FROM DEVICE PACKAGING

|                                                                                       |                                                                                                                                                                                                                                                                                                             |
|---------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|    | MEDICAL DEVICE                                                                                                                                                                                                                                                                                              |
|    | CONTAINS BIOLOGICAL MATERIAL OF ANIMAL ORIGIN                                                                                                                                                                                                                                                               |
|    | KEEP AWAY FROM SUNLIGHT                                                                                                                                                                                                                                                                                     |
|    | KEEP DRY                                                                                                                                                                                                                                                                                                    |
|    | UPPER LIMIT OF TEMPERATURE, 25°C, 75°F, LOWER LIMIT OF TEMPERATURE, 32°F, 0°C                                                                                                                                                                                                                               |
|    | CONSULT INSTRUCTION FOR USE                                                                                                                                                                                                                                                                                 |
|    | DO NOT USE IF PACKAGING IS DAMAGED                                                                                                                                                                                                                                                                          |
|    | DO NOT REUSE                                                                                                                                                                                                                                                                                                |
|   | CAUTION                                                                                                                                                                                                                                                                                                     |
|  | STERILIZED USING ETHYLENE OXIDE                                                                                                                                                                                                                                                                             |
|  | DISTRIBUTOR                                                                                                                                                                                                                                                                                                 |
|  | DO NOT RE-STERILIZE                                                                                                                                                                                                                                                                                         |
|  | SINGLE STERILE BARRIER SYSTEM                                                                                                                                                                                                                                                                               |
|  | DOUBLE STERILE BARRIER SYSTEM                                                                                                                                                                                                                                                                               |
|  | MANUFACTURER                                                                                                                                                                                                                                                                                                |
|  | CATALOG NUMBER                                                                                                                                                                                                                                                                                              |
|  | UNIQUE DEVICE IDENTIFIER                                                                                                                                                                                                                                                                                    |
|  | LOT NUMBER                                                                                                                                                                                                                                                                                                  |
|  | EXPIRY DATE                                                                                                                                                                                                                                                                                                 |
|  | CE MARK, NOTIFIED BODY: DEKRA CERTIFICATION B.V.                                                                                                                                                                                                                                                            |
|  | KERECIS LIMITED<br>Eyrargata 2<br>400 Isafjördur<br>Iceland<br>Phone: +354 419 8000<br>Email: <a href="mailto:info@kerecis.com">info@kerecis.com</a><br><a href="http://www.kerecis.com">www.kerecis.com</a><br><br>Patented: US 8,613,957.<br>Kerecis and the Kerecis logo are trademarks of Coloplast A/S |
|  | Notified Body: DEKRA Certification B.V.                                                                                                                                                                                                                                                                     |

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